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Investigation of analysis for vulcanized natural rubber by
inverse measurements with FG-MAS probe
(Application & Research Group, JEOL Ltd.¹,
Department of Chemistry, Nagaoka University of Technology²)
Satoshi Sakurai¹, Hiroaki Utsumi¹, Seiichi Kawahara²

Analysis of crosslinking junctions of crosslinked rubber samples has been made through solid-, solution- and latex-state NMR spectroscopies up to date. If we apply ordinary pulse techniques of solution-state NMR spectroscopy under MAS, it is possible to positively analyze the crosslinking junctions, more exactly, because the rubber samples have high mobility. In this case, the analysis through inverse measurements are quite effective to shorten a time for the measurement. However, because the signals assigned to the main chains are too big, the huge correlation signals of the main peaks interfere with the small signals of the crosslinking junctions under a resolution of the ¹H observation. In the present study, an attempt to analyze the crosslinking junctions of the rubber samples was made by solution-state NMR spectroscopy under MAS with selective neglect of the main correlation signals.

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
Fig. 1 Pulse Sequence of yselect ¹H-¹³C HMQC

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